

ALBAcyte®REAGENT RED BLOOD CELLS FOR IDENTIFICATION OF UNEXPECTED ANTIBODIES
Antibody Identification (16-Cell) - [REF] Z473U



Patient/Sample Details				
Name:	D.O.B:	Sample Identification:	ABO Rh:	
Diagnosis History:	Hospital:	DAT:	Anti-IgG:	Anti-C3:

Lot No: V225877
Expiry Date: 2020.11.09 U.S License 1807

Cell #	Rh-hr	Donor	Rh-hr							Kell						Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens			TEST RESULTS				Cell #		
			D	C	E	c	e	f	*V	C ^w	K	k	Kp ^a	Kp ^b	*Js ^a	*Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	M	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wr ^a	Special Types						
1	R ₁ ^w R ₁	2768020435687	+	+	0	0	+	0	NT	+	0	+	0	+	0	+	+	+	0	+	0	+	+	0	+	0	+	0	+	0	+	*0	0					1
2	R ₁ R ₁	DN00390508	+	+	0	0	+	0	NT	0	+	+	0	+	0	+	0	+	+	0	+	0	+	+	0	+	0	+	+	0	+	+	0					2
3	R ₂ R ₂	6106011312018	+	0	+	+	0	0	NT	0	0	+	0	+	0	+	+	+	+	0	+	+	+	+	+	0	+	0	+	0	0					3		
4	r'r	1095603	0	+	0	+	+	+	NT	0	0	+	0	+	0	+	+	0	+	0	+	0	+	0	+	0	+	0	+	0	0					4		
5	r'r	4363184	0	0	+	+	+	+	NT	0	0	+	0	+	0	+	+	+	+	0	+	+	+	0	+	+	+	+	+	+	0					5		
6	R ₀ r	1371316	+	0	0	+	+	+	NT	0	0	+	0	+	+	+	0	0	+	0	0	+	+	0	+	+	0	+	0	0					6			
7	rr	4389710	0	0	0	+	+	+	NT	0	+	+	0	+	0	+	+	0	+	0	+	+	0	+	0	+	0	+	0	0					7			
8	rr	174097	0	0	0	+	+	+	NT	0	+	+	0	+	0	+	0	+	0	+	+	+	0	+	+	0	+	0	+	*	+	0					8	
9	rr	1300117	0	0	0	+	+	+	NT	0	0	+	0	+	0	+	0	+	+	0	+	+	0	+	+	0	+	0	+	+	0					9		
10	rr	3336181	0	0	0	+	+	+	0	0	0	+	+	0	+	0	+	+	0	+	0	+	0	+	0	+	0	+	+	0					10			
11	R ₁ R ₂	635243414	+	+	+	0	+	0	NT	0	0	+	0	+	0	+	0	+	+	0	+	+	0	+	0	+	0	+	+	0					11			
12	R ₁ R ₁	890920	+	+	0	0	+	0	NT	0	0	+	0	+	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	0					12			
13	R ₁ r	523326	+	+	0	+	+	+	NT	0	+	0	0	+	0	+	+	+	0	0	+	0	+	+	+	+	+	+	+	0					13			
14	R ₂ R ₂	77430	+	0	+	+	0	0	NT	0	+	+	0	+	0	+	+	+	0	+	0	+	0	+	0	+	0	+	+	0					14			
15	R ₂ R ₂	53162855	+	+	+	+	0	0	NT	0	0	+	0	+	0	+	+	+	0	+	+	+	+	0	+	0	+	*	0	0					15			
16	rr	11073434	0	0	0	+	+	+	NT	0	0	+	0	+	0	+	+	+	0	0	+	+	+	0	S+	0	+	+	0					16				
Patients cells																																						

Notes:
1. All cells are DAT negative.
2. 's' = strong reactor, 'w' = weak reactor, NT = not tested
3. The f antigen status has been determined presumptively
4. Cell #6 is designated as R₀r but the phenotyping test results cannot exclude the possibility of the red blood cell genotype being R₀R₀.
* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.
Quotient, 301 South State Street, S-204, Newtown, PA 18940, USA
Product Technical Support Tel: 1-888-228-1990 Web: www.quotientbd.com

Conclusions/Further Testing Required:	
Signature/Initials:	Date: