

ALBAcyte® REAGENT RED BLOOD CELLS FOR USE IN IDENTIFICATION OF UNEXPECTED ANTIBODIES

REF Z471U & Z472U



Patient/Sample Details				
Name:	D.O.B:	Sample Identification:	ABO Rh:	
Diagnosis History:	Hospital:	DAT:	Anti-IgG:	Anti-C3:

Lot No:	V248675 & V248678
Expiry Date:	2022.07.18
U.S License	1807

Cell #	Rh-hr	Donor	Rh-hr								Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens		TEST RESULTS				Cell #		
			D	C	E	c	e	f	*V	C ^w	K	k	Kp ^a	Kp ^b	*Js ^a	*Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	M	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	W ^r ^a	Special Types				
1	R ₁ ^w R ₁	4141621	+	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	+	0	+	+	+	+	+	0	+	+	+					1
2	R ₁ R ₁	52005658	+	+	0	0	+	0	NT	0	+	+	0	+	0	+	+	0	+	+	0	+	+	+	0	+	+	0	+	*+	0					2
3	R ₂ R ₂	1352908	+	0	+	+	0	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	0	+	+	0					3	
4	r'r	10839983	0	+	0	+	+	+	0	0	0	+	0	+	0	+	+	+	+	0	+	+	+	+	0	0	0	+	+	0					4	
5	r''r	3097391	0	0	+	+	+	+	NT	0	0	+	0	+	0	+	+	+	+	0	+	+	+	+	0	0	0	+	0	0					5	
6	R ₀ r	54742364	+	0	0	+	+	+	NT	0	0	+	0	+	0	+	0	0	+	+	0	0	+	+	0	+	+	0	+	*+	0					6
7	rr	54214044	0	0	0	+	+	+	NT	0	+	+	0	+	0	+	+	+	+	0	0	+	+	+	0	+	0	+	*+	*0					7	
8	rr	1580417	0	0	0	+	+	+	NT	0	0	+	+	+	0	+	+	+	+	0	+	+	+	+	+	0	+	+	+	0					8	
9	rr	4256454	0	0	0	+	+	+	0	0	0	+	0	+	0	+	+	0	+	+	0	+	+	+	+	+	+	+	0	0					9	
10	rr	3146843	0	0	0	+	+	+	0	0	+	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	0					10	
Patient Cells																																				

Notes:

- 1. All cells are DAT negative.
 - 2. M, N, S, s, Fy^a and Fy^b antigens are destroyed or depressed following enzyme treatment.
 - 3. 's' = strong reactor, 'w' = weak reactor, NT = not tested
 - 4. The f antigen status may have been determined presumptively based on the Rh-hr phenotype.
- * Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Conclusions/Further Testing Required:
Signature/Initials: _____ Date: _____