

Expected Test Results	Reagents					Reagent Red Cells						QC Kit lot Number:				
	Anti-A	Anti-B	Anti-AB	Anti-D	Rh Control	A1	A2	B	O	Screen Cell 1		Screen Cell 2		Screen Cell 3		Expiration Date:
										37°	IAT	37°	IAT	37°	IAT	
QC Cell 1 – AB rr	1-4+	1-4+	1-4+	0	0											
QC Cell 2 – O R1r	0	0	0	1-4+	0											
QC Antibody						1-4+	1-4+	1-4+	0	0-2+	1-4+	0-2+	1-4+	0-2+	1-4+	

Daily Test	Reagents					Reagent Red Cells						Additional Reagents						
	Anti-A	Anti-B	Anti-AB	Anti-D	Rh Control	A1	A2	B	O	Screen Cell 1		Screen Cell 2		Screen Cell 3		AHG	Enhancement Media	IgG Sensitized Cells
										37°	IAT	37°	IAT	37°	IAT			
Lot No:																Lot #:	Lot #:	Lot #:
Expiry Date:																Exp:	Exp:	Exp:
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																Tech:	Date:	
Unsatisfactory*																		

Lot No:																Lot #:	Lot #:	Lot #:
Expiry Date:																Exp:	Exp:	Exp:
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																Tech:	Date:	
Unsatisfactory*																		

Lot No:																Lot #:	Lot #:	Lot #:
Expiry Date:																Exp:	Exp:	Exp:
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																Tech:	Date:	
Unsatisfactory*																		

Lot No:																Lot #:	Lot #:	Lot #:
Expiry Date:																Exp:	Exp:	Exp:
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																Tech:	Date:	
Unsatisfactory*																		